

## IV INFUSION BY CURLIN PUMP INSTRUCTIONS



602-997-0500

Gather all the supplies you will be using:

- a. Alcohol Wipes
- b. Medicine Bag
- c. IV Tubing
- d. Saline Flush

**Wash your hands with antibacterial soap and water.**

### ***Installing Batteries:***

1. Locate battery door on back of pump.
2. Open the lock; slide door open in downward motion.
3. Install two "C" batteries into pump and close door by sliding door back into place.

### ***Loading IV tubing:***

1. Fully open the door located on the top of the pump.
2. Remove IV tubing from package and locate the **Blue** Tubing-Guide and **Yellow** Flow stop.
3. Insert the **Blue** Tubing-Guide into the receptacle on the right side of the pump (look for the blue arrow).
4. Twist tab off **Yellow** Flow-Stop.
5. Insert **Yellow** Flow-Stop into the receptacle on the left side of the pump (look for the yellow arrow).
6. Close pump door and latch it securely.
7. Remove protective cap from spiked end of tubing.
8. Push spiked end of tubing into IV bag using a twisting/pushing motion.

### ***Starting Infusion for the first time:***

1. Turn pump on (locate **ON/OFF** button on the top left of keypad);
  - a. Pump will sound 3 beeps and go through a series of self tests on the display screen.
  - b. When pump has completed its self tests the pump then will display the following:

ENTER TO SELECT  
PROGRAM  
BIOMED SETUP

2. Press the yellow **YES** button located on pump keypad to select “PROGRAM”.

3. The pump will then display the following:

ENTER TO SELECT  
*(RESUME)\*\**  
REPEAT RX  
NEW PROGRAM

4. Press the yellow **YES** button located on the pump keypad to select “REPEAT RX”.

- a. **\*\***This option will show up the next time you will change the IV bag. Follow the same steps as listed. **DO NOT** choose this option. Always choose “REPEAT RX” after you have attached a new IV bag. Only choose “RESUME” if you have **paused** the infusion for whatever reason.

5. The pump will then run through its pre-programmed settings and will then display:

RUN TO START  
NO TO REVIEW  
OPTIONS OR PRIME

4. Press and release the **PRIME/BOLUS** button located at the lower right corner of pump keypad.
5. Press and HOLD the **PRIME/BOLUS** button until tubing is completely filled with fluid.
6. When priming is complete press **YES** button to exit PRIME screen.
7. Wipe tip of line with alcohol prep pad (approx. 10-15 seconds).
8. Flush IV catheter using **Saline**.
9. Connect IV tubing to tip and then Push the **RUN/PAUSE** button to start infusion.